


## INCREMENTAL ANIMAL HANDLING SYSTEM



A complete integrated protection system for animal handlers and their subjects. To safeguard personnel while handling animals, Lomir manufactures a complete range of handling equipment designed and fabricated for flexibility, comfort and safety during long periods of wear. No other company offers such a range of sizes; Small Medium, Large and X-Large. Available as pairs or singles. For more information, or to order please visit us online at [www.lomir.com](http://www.lomir.com) or call 1-877-425-3604.



### LEVEL 1

Kevlar® Gloves - 100% fine-knit Kevlar® gloves. May be worn on their own for handling small animals, or under handling gloves for species such as primates.



### LEVEL 2

Kevlar® Arm Protectors - Arm protectors made from 100% fine-knit Kevlar® with sewn-in reinforcement. Designed to be worn either on their own or beneath a pair of animal handling gloves.



### LEVEL 3

Leather Gauntlets - Full fingered, full length gauntlet. Glove section made from 100% cowhide, elbow length cuffs made from split grain cowhide. Ruggedly constructed with internal seams throughout for optimum protection and durability.



### LEVEL 4

Leather Doublers - For the highest level of protection, a 100% cowhide doubler is available with heavy duty snaps for attachment over the handling gauntlets. Full thumb and index finger closure.



458 East Main Street, Malone, NY, 12953  
Tel: (518) 483-7697 Fax: (518) 483-8195  
Toll Free: (877) 425-3604  
[www.lomir.com](http://www.lomir.com)

95 Huot, Notre-Dame-de-l'Île-Perrot, Qc., J7V 7M4  
Tel: (514) 425-3604 Fax: (514) 425-3605  
Toll Free: (877) 425-3604  
[www.lomir.com](http://www.lomir.com)

