

## DERMAL APPLICATIONS INSERTS AND JACKETS



Designed for studies where test articles are applied by the dermal route, this system provides a non-occlusive covering for the application site, reducing the possibility of oral ingestion, while allowing in-life investigations to be carried out with ease. Used in regulatory toxicology studies, these dermal inserts and jackets are suitable for many products that are administered by this route.

For more information, or to order please visit us online at [www.lomir.com](http://www.lomir.com) or call 1-877-425-3604.



Stock and custom dermal jackets are available for all subjects studied in dermal procedures.

Lomir's animal jackets and inserts are manufactured for a range of species such as rabbits, guinea pigs and rodents, with custom models available for dogs and swine.



Quality materials and construction ensure durability and efficiency over repeated usage. Adjustable velcro connecting strip provides proper fit and ease of use.



458 East Main Street, Malone, NY, 12953  
Tel: (518) 483-7697 Fax: (518) 483-8195  
Toll Free: (877) 425-3604  
[www.lomir.com](http://www.lomir.com)

95 Huot, Notre-Dame-de-l'Île-Perrot, Qc., J7V 7M4  
Tel: (514) 425-3604 Fax: (514) 425-3605  
Toll Free: (877) 425-3604  
[www.lomir.com](http://www.lomir.com)

