

PRIMATE CAPTURE NETS



Lomir's capture nets are designed to restrain animals while protecting both the subject and the handler. They are suitable for both new and old world primates. The net and pole are available in multiple combinations to match the size, weight and strength of the subject. Replacement nets are also available. As with all Lomir's products, custom specifications are also undertaken.

For more information, or to order please visit us online at www.lomir.com or call 1-877-425-3604.



The net is a robust, fray-resistant polyester mesh fabric, double sewn and lock stitched for strength and durability. A reinforced channel is sewn in to accommodate a thick braided cord. This can be tightened once the subject has been netted for additional constraint and safety.



Capture nets are suitable for both new and old world primates.



Two energy absorbing grips are ergonomically positioned on the handle for security and comfort.



Corners of the stainless steel hoop are rounded for additional safety.



458 East Main Street, Malone, NY, 12953
Tel: (518) 483-7697 Fax: (518) 483-8195
Toll Free: (877) 425-3604
www.lomir.com

95 Huot, Notre-Dame-de-l'Île-Perrot, Qc., J7V 7M4
Tel: (514) 425-3604 Fax: (514) 425-3605
Toll Free: (877) 425-3604
www.lomir.com