

## SMALL ANIMAL INFUSION SYSTEM



Completely integrated rodent infusion systems providing maximum procedural flexibility and precision. Infusion sets may be purchased complete or as individual components and are also available pre-serilised. For more information, or to order please visit us online at [www.lomir.com](http://www.lomir.com) or call 1-877-425-3604.



### SWIVELS

High impact thermoplastic, precision stainless steel hypodermic tubing, and Teflon seals. Single or double fluid channels in a range of needle sizes.



### JACKETS

Made of double layers of Lycra/Spandex fabric with Velcro fasteners for adjustment and secure fit. Available in 5 standard sizes, (X-small, Small, Medium, Large and X-large), order by weight range of your animals.




### TETHERS

Manufactured from a strong but light stainless steel spring, in three standard lengths: 12" (30cm), 18" (46cm) and 24" (61cm) for single or multiple channel use. The standard anchor is a molded button with a sewn Velcro tab for adjustable fit when using the jackets.



### TAIL CUFF

Constructed of stainless steel with self-threading tether attachment point. Used as an alternative to the Rodent Jacket.




### CONNECTORS

These provide a secure attachment between Swivels and Tethers ensuring an ease of movement through 360°.



### RETAINERS

Stainless steel and high-impact plastic construction. Available in two designs selected by the type of cage in use.



458 East Main Street, Malone, NY, 12953  
Tel: (518) 483-7697 Fax: (518) 483-8195  
Toll Free: (877) 425-3604  
[www.lomir.com](http://www.lomir.com)

95 Huot, Notre-Dame-de-l'Île-Perrot, Qc., J7V 7M4  
Tel: (514) 425-3604 Fax: (514) 425-3605  
Toll Free: (877) 425-3604  
[www.lomir.com](http://www.lomir.com)