

SMALL PRIMATE HANDLING

When handling small animals such as New World primates dexterity coupled with protection is vital for the security of both the subject and handler. Lomirs' small primate handling gloves offer a smart solution for protection during long periods of wear.

Manufactured from a soft yet durable deerskin, these gloves are adjustable for a secure and comfortable fit and lined with knitted Kevlar to provide extra protection and slash resistance.



Kevlar lining provides extra protection from slashing and tearing.



Adjustable wrist strap allows for an accurate comfortable and secure fit.



SAFE, DURABLE AND COMFORTABLE

Sizes Available: Small, Medium, Large, X-Large

Cat #	Description
HG 011KS	Deer skin handling gloves with Kevlar lining. Small
HG 011KM	Deer skin handling gloves with Kevlar lining. Medium
HG 011KL	Deer skin handling gloves with Kevlar lining. Large
HG 011KX	Deer skin handling gloves with Kevlar lining. X-large



1. Deerskin leather, thin yet strong and durable
2. Lined with Kevlar for slash resistance
3. Adjustable Velcro wrist strap
4. Soft and form fitting will not crack or harden after wetting



HANDLING MADE EASY



Lomir
BIOMEDICAL INC.

99 East Main Street, Malone, NY, 12953
Tel: (518) 483-7697 Fax: (518) 483-8195
Toll Free: (877) 425-3604
www.lomir.com

95 Huot, Notre-Dame-de-l'Île-Perrot, Qc., J7V 7M4
Tel: (514) 425-3604 Fax: (514) 425-3605
Toll Free: (877) 425-3604
www.lomir.com