

RODENT SLING

SLING FRAME & COVER




Figure 1: Rodent Sling Frame & Cover

FRAME

The rodent sling frame is routinely supplied with a one-size adjustable canvas and mesh cover. This is shown in the picture above (Figure 1). The design allows for animals to be quickly located in the sling and for easy adjustment and fit between individual animals.

SLING SUIT

The rodent sling suit (purchased separately) is best described as a full body rat jacket with canvas extensions through which the supports of the frame are passed. The body portion of the suit is made from our lycra/spandex rat jacket fabric and closes by means of a Velcro insert. This provides for a quick, secure fit and adjustment between individual animals. The extensions are made from canvas fabric with a sleeve on each side through which the bars of the frame are passed. The overall dimensions make the suit compatible with the sling frame and the jacket portion of the suit corresponds to a full body rat jacket.



Should you have any questions or need more information please contact our Customer Service Department toll free in North America at 1-877-425-3604 or send an e-mail to info@lomir.com.

Lomir Biomedical Inc.

95 Huot, Notre-Dame-de-l'Île-Perrot, QC
Canada J7W 1M3
Tel: 514-425-3604
Toll Free: 1-877-425-3604
Fax: 514-425-3605
Email: info@lomir.com

Lomir Biomedical Inc.

C/O Agenda Life Sciences
UK PO Box 24 Hull
HU128YK, UK
Tel: 03456 44 55 45
Email: info@lomir.com

Lomir Inc.

PO Box 778
Malone, NY
USA 12953
Tel: 518-483-7697
Toll Free: 1-877-425-3604
Email: info@lomir.com